NorthEast Cryo/Bio/Science
Home of CryoTooth & Store-A-Tooth

CryoTooth Collection Kit Instructions

Your kit has everything needed to safely send baby or wisdom teeth to our lab. Follow these steps, or give this sheet to your dentist or oral surgeon.

When your kit arrives

Freeze the cold packs right away
Take them out of the kit and put them in your freezer. They need to be frozen for at least 24 hours before the tooth is extracted.

Save the insulated return shipper
The insulated cold-shipper box your kit arrives in is what you'll use to send the tooth back.

A return label is included
You'll attach it to the insulated return shipper when you send the tooth back.

If a dental professional is collecting the tooth: bring all kit components to the appointment, including the documents and the frozen cold packs, and ask them to follow the steps below. Plan to drop the kit at a staffed UPS location before its latest cutoff for next-day delivery.

When the tooth comes out



1 **Place the tooth or teeth in the Save-A-Tooth-style container** and secure the cap.



2 **Place the container in the protective specimen bag** and close it tightly. *(The bag is marked "biohazard.")*




3 **Place the sealed bag at the bottom of the styrofoam cooler.**




4 **Surround the container with the frozen cold packs.**



5 **Place the cooler in the insulated return shipper and seal the cooler with its cover.** Don't forget the Donor Info Card.



6 **Seal the insulated return shipper with packing tape** and attach the return label.



7 **Bring the kit to a staffed UPS location the same day**, before the overnight shipping deadline.

Shipped? Let us know.

Email info@NECryoBioScience.com so we can watch for its arrival. Questions? Call 888-441-GUMS (4867).